EN

Cryogenic Vials are designed for storing biological materials for in vitro diagnostic. Single use device used by Laboratory professionals only. Sterilized by gamma radiation at a SAL of 10^{-3} . They can withstand temperatures from +121°C down to -196°C.

INSTRUCTIONS FOR USE

1. Identify sample with appropriate laboratory marking pen.
2. Fill the Cryovial® with a maximum of 90% of the tube's volume to permit sample expansion.
3. Verify that the Cryovial® and screw cap are completely dry before closing; liquid drops will impair the seal in the gas phase of liquid nitrogen.
4. Tighten the Cryovial® screw cap normally;
5. Store the vials in the gas phase of liquid nitrogen.
6. To thaw, place the Cryovial® in a water bath (37°C to 40°C), keeping the vials in constant motion. Remove the Cryovial® from the water bath when the ice has melted.

Warning:

- Graduation markings on tubes are for approximate volume and are not a calibrated measuring instrument.
- Over-tightening the cap may squeeze the silicone washer between the tube and cap and could increase the risk of leakage.
- Do not use if package is damaged and sterility is compromised.
- **Use only in the gas phase of Liquid Nitrogen. Immersion in the liquid phase may cause infiltration into the vials and could result in leaking or shattering when brought to room temperature.**
- Use safety equipment: gloves, hoods, face shields, etc. when removing vials from liquid nitrogen containers.
- Safely dispose of infectious substances following applicable rules and regulations on waste disposals.
- Please report any serious incident to Simport and to the competent authority of the Member State where you are established.

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