

Product Name: Specimen Container

Catalogue No.: C581-300

This document replaces any previous version

1. Product Description:

- **300 ml Specimen Container:** Disposable and non-sterile with tamper evident lid

2. Packaging:

- **Case:** 100 units per case

3. Product Specifications:

- **Material:**
 - Container: Polypropylene
 - Lid: Polypropylene
- **Temperature range:** -20°C to +80°C

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
- **CONEG / RoHS :** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals

5. Quality Assurance:

- No presence of contamination in plastic
- Visual attributes
- Closure verification

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging

8. Recommended Use:

- Check cap closure when using biohazard material and / or chemical reagents
- For freezing, check volume to prevent overflowing

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*