

Product Name: Disposable Funnel

Catalogue No.: F490-1

This document replaces any previous version

1. Product Description:

- **Funnel:** Disposable and non-sterile

2. Packaging:

- **Case:** 4 stacks of 25 units / 100 units per case

3. Product Specifications:

- **Material:** Polystyrene
- **Top Inner Diameter:** 57 mm
- **Total Length:** 69.5 mm
- **Use with Paper Diameter:** 11 cm
- **Tapered at a 60° angle with inside fluting**

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** plastics conform to FDA 21 CFR 177-1640
- **CONEG:** Plastics and colorants are in conformity with CONEG standards for heavy metals

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Careful for spills when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*