

**Product Name: Biopsy Foam Pads**

**Catalogue No.: M476-4**

**This document replaces any previous version**

## 1. Product Description:

- **Biopsy Foam Pad:** Non-sterile, disposable rectangular foam pad for the Simport Slimsette cassettes

## 2. Packaging:

- **Case:** 10 bags of 1 000 units / 10 000 units per case

## 3. Product Specifications:

- **Material:**
  - Blue Polyurethane Foam
- **Dimensions:** 27.4 x 25.4 x 2 mm
- **Temperature range:** -40°C to +121°C
- **Autoclavable at +121°C for up to 30 minutes**

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

## 5. Quality Assurance:

- No presence of contamination in foam
- Visual attributes
- Dimension measurements

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

## 8. Recommended Use:

- Fits Simport Slimsette cassettes
- For use to hold biopsies inside histology cassettes
- check cassette's lid closure before tissue processing

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*