

**Product Name:** Write-on Marker pen

**Catalogue No.:** M495-12BK

This document replaces any previous version

**1. Product Description:**

- **Marker Pen:** Laboratory pen for marking on embedding cassettes and slides, Non-sterile,

**2. Packaging:**

- **Case:** 12 units per case

**3. Product:**

- **Material:** Permanent, black ink marker pen

**4. Standards and Conformity:**

- **ISO 2859-1:** Sampling and inspection procedures

**5. Quality Assurance:**

- Resistant to most histology solvents such as Xylene, Alcohol, Acetone and Formalin and xylene substitutes.
- Writing will remain clean, clear and legible throughout the range of laboratory procedures, including processing and staining.
- Ink will dry quickly and will not run or smear
- Ink will last through fixatives, baths, dips and stains

**6. Traceability:**

- **Lot No. Composition:** 8 or 9 digits  
**The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

**7. Storage Conditions:**

- Store at room temperature in normal warehouse conditions
- Store in upright position

**8. Recommended Use:**

- Used for writing information on cassettes and slides

**OWNERSHIP OF MATERIALS:** Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.