

Product Name: **DissecTable™**

Catalogue No.: **M620**

This document replaces any previous version

1. Product Description:

- **M620 Dissecting board:** Printed on both sides. Large and small cutting areas including dimensional scales in inches and centimeters
- **Compatible with:**
 - **M625** Board base (sold separately) made of Polystyrene, intended to hold the dissecting board and will retain excess fluids when necessary.

2. Packaging:

- **Case:** 1 individually wrapped unit per case

3. Product Specifications:

- **Material:**
 - High-Density Polyethylene
- **Temperature range:**
 - Base -80° C TO + 121° C
- **Dimensions:** 23 x 16 x 0.5 in.

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging or product

OWNERSHIP OF MATERIALS:

Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.