

Product Name: Dispocut™

Catalogue No.: M630

This document replaces any previous version

1. Product Description:

- **M630 Dissecting board:** Printed on both sides. Large and small cutting areas including dimensional scales in inches, centimeters and millimeters. Offered in three sizes

2. Packaging:

- M630-1 4 packs of 24 units / 96 per case
- M630-2 4 packs of 12 units / 48 per case
- M630-3 2 packs of 12 units / 24 per case

3. Product Specifications:

- Material:
 - Polypropylene
- Dimensions:
 - M630-1: 152 x 203 mm
 - M630-2: 229 x 305 mm
 - M630-3: 305 X 483 mm
- Temperature range: -27°C to +71°C

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging or product

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*