

Product Name: ChillBlock®

Catalogue No.: S700-35

This document replaces any previous version

1. Product Description:

- **Thermal Conductive Rack:** Non- Sterile rack has 9 wells for 15 ml centrifuge tubes

2. Packaging:

- **Case:** case of 1 unit

3. Product Specifications:

- **Material:** Metal alloy
- **Fits 15 ml (Simport T408) to 15 ml Centrifuge Tubes**
- **Outer dimensions:** 3.1 x 3.1 x 4.2 in (7.9 x 7.9 x 10.7 cm)

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **Thermo-conductive:** rapidly adapts to temperatures from -150°C to +100°C

5. Quality Assurance:

- Visual attributes

6. Traceability:

- **Lot No. Composition:** 8 to 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label;
 2. On label inserted inside the master case;
 3. On the inner bag.

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging

8. Recommended Use:

- None

OWNERSHIP OF MATERIALS:

Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.