

Product Name: **Deep Well Plates**

Catalogue No.: **T110-25**

This document replaces any previous version

1. Product Description:

- **Matt Cover for 1.2ml 96 Deep Well Plates:** Non-sterile, designed to fit the *Simport* Bioblock™ T110-5

2. Packaging:

- **Case:** 24 units per sealed bag

3. Product Specifications:

- **Material:**
 - Matt Cover : Polyolefin Elastomer
- **Temperature range:** -80°C to +80°C

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension and volume measurements.

6. Traceability:

- **Lot No. Composition :** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any damage to the packaging

8. Recommended Use:

- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*