



M520 and M521-SL

Microscreen™ I

Small Biopsy processing/embedding cassettes
The finest design for small and delicate specimens

The Microscreen™ cassette is made of a strong nylon mesh and designed to handle small, fragile biopsies.

- Sturdy nylon mesh improves fluid exchange by reducing surface tension of reagents
- Accelerates fixation/preparation of tissue samples
- Eliminates the need for pads, sponges, nylon or lens paper bags
- Reduces number of incidents involving non-fixed or insufficiently fixed tissues
- Inner chamber design for greater biopsy protection
- Large front tab for easier lid opening
- The hinge design configuration allows for multiple opening and closing



Specifications

- Cassette frame made of acetal, resistant to histological solvents
- Mesh made of nylon: 1500 square openings (300 µm each)
- Eight 11 mm slots in base and lid to help cassette sink rapidly
- Anterior writing area at 45° angle

Bulk			In QuickLoad™ sleeve		
Cat. No.	Qty/Bx	Qty/Cs	Color	Cat. No.	Qty/Cs
M520-2	250	1000	○ White	M521-2SL	750
M520-3	250	1000	● Pink	M521-3SL	750
M520-4	250	1000	● Green	M521-4SL	750
M520-5	250	1000	● Yellow	M521-5SL	750
M520-6	250	1000	● Blue	M521-6SL	750
M520-7	250	1000	● Peach	M521-7SL	750
M520-8	250	1000	● Tan	M521-8SL	750
M520-9	250	1000	● Gray	M521-9SL	750
M520-10	250	1000	● Lilac	M521-10SL	750
M520-11	250	1000	● Orange	M521-11SL	750
M520-12	250	1000	● Aqua	M521-12SL	750

