

Product Name: Caps for Transport Tube

Catalogue No.: T552WTP

This document replaces any previous version

1. Product Description:

- **Tamper Evident Cap:** Disposable, non-sterile lip seal cap for 5ml, 10ml and 12ml transport tubes.
- **Compatible with:**
 - T552-5ATTP 5ml Transport Tube Series
 - T550-10ATTP 10ml Transport Tube Series
 - T550-12ATTP 12ml Transport Tube Series

2. Packaging:

- **Case:** 1 000 units per pack

3. Product Specifications:

- **Material:**
 - High Density Polyethylene
- **Color:** White
- **Temperature range:** -90°C to +100°C

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG/RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free
- **BSE / TSE:** Material is BSE / TSE Free

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimensional measurements
- Closure verification
- Leak proof testing in vacuum at 71.3cm Hg (95 kPa tested)

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*