

Product Name: **Deep Base Mold**

Catalogue No.: **M475-10**

This document replaces any previous version

1. Product Description:

- **Base Mold:** Non-sterile, disposable deep base mold, designed for M512 MACROSETTE™

2. Packaging:

- **Case:** 500 units per case

3. Product Specifications:

- **Material:**
 - PVC
- **Temperature range:** -30°C to +50°C
- **Autoclavable:** No

Compartment size: 37mm X 24mm X 10mm
Volume: 5 ml

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

5. Quality Assurance:

- No presence of contamination on stainless steel
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Accessory to be used in combination with cassettes for embedding procedures

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*