

Product Name: UriFlex®

Catalogue No.: B351-4

This document replaces any previous version

1. Product Description:

- **4 Liters, 24-Hr Urine Container:** Disposable, Non-sterile assembled with screw cap with pour spout

2. Packaging:

- **Case:** 30 units per case

3. Product Specifications:

- **Material:**
 - **Bottle: laminated:**
 - Outer film PET
 - Middle film Nylon
 - Inner liner PE
 - **Cap: Polypropylene**
- **Temperature range: 2 to 8°C**

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

5. Quality Assurance:

- No presence of contamination in plastic;
- Visual attributes;
- Closure and valve verification.

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case;
 3. On the inner bag.

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*