

Product Name: Deep well plate

Catalogue No.: T110-49

Material:

Natural Polypropylene *inquire for colors*

Description:

- Deep well plate with 384 V bottom square wells
- Volume 240 uL
- Dimensions: SBS - American Standards (ANSI) (ANSI/SBS 1-2004), see drawing for specifics

Packaging:

Qte/pk: 10 units

Qte/Case: 50 units

Standards and Conformity:

Temperature range: -86°C to +121°C

Autoclavable: +121°C for up to 20 minutes. *The dimensional stability of the disposables can be affected.*

Centrifuge: resistance up to 4000 g

RNase, DNase, DNA and Pyrogen free: certified

ISO 2859-1: Sampling and inspection

FDA: Plastics conform to FDA 21 CFR 177.1520

CONEG / RoHS: Plastics comply

REACH: Plastic comply

LATEX: Plastic is Latex free

BSE/TSE: Plastic is BSE / TSE Free

Produced and packed in an ISO class 8 clean environment

Recommended Use: In automated equipment please follow the equipment manufacturer's instructions

Storage Conditions: Store at room temperature in normal warehouse conditions. Protect from any possible contamination, direct sunlight and UV light, and any damage to the packaging.

Disclaimer: *The information contained in this document shall not constitute a representation or warranty of the suitability for any specific purpose. To the best of our knowledge, the information contained is reliable and true based on our internal and external testing as well as statements supplied to us by third party suppliers.*