

Product Name: The Multirack™

Catalogue No.: S600-30

1. Product Description:

- **Rack:** Non-Sterile, 18 position rack
- **Available Colors:**
 - Blue, Yellow, Green, Orange, Lilac

2. Packaging:

- **Case:** 10 units per case

3. Product Specifications:

- **Material:**
 - Acetal
 - Not Autoclavable
 - Total capacity: Holds 18 tubes of 25 to 30 mm

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.2470
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions

8. Recommended Use:

Rack for holding tubes

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*