

**Product Name:** The Multirack™ Jr.

**Catalogue No.:** S610-16

## 1. Product Description:

- **Rack:** Non-Sterile, 30 position rack
- **Available Colors:**
  - Blue, Yellow, Green, Orange, Lilac

## 2. Packaging:

- **Case:** 10 units per case

## 3. Product Specifications:

- **Material:**
  - Acetal
  - Not Autoclavable
  - Total capacity: Holds 30 tubes up to 16 mm

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.2470
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
  
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions

## 8. Recommended Use:

Rack for holding tubes

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*