

**Product Name:** **Specimen Container**

**Catalogue No.:** **C570-12**

This document replaces any previous version

**1. Product Description:**

- **4 oz Specimen container:** Disposable, Non-sterile Specimen container with separate lids.

**2. Packaging:**

- **Case:** 5 packs of 100 units / 500 units per case

**3. Product Specifications:**

- **Material:**
  - Tube: Polypropylene
  - Cap: Polyethylene
- **Temperature range:** -90°C to +100°C
- **Tubes have molded-in graduations up to 128 ml**

**4. Quality Assurance:**

- Clear, no presence of contamination in plastic
- Visual attributes
- Closure verification

**5. Traceability:**

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

**6. Storage Conditions:**

- Store at room temperature in normal warehouse conditions

**7. Recommended Use:**

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations

**OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*