

Product Name: **CytoSep™**

Catalogue No.: **M968-1**

1. Product Description:

- **Cytology funnel :** Disposable cell concentrator for 50 to 400ul sample volume. Filters are sold separately. Cell concentration area: 7mm diameter
- **Compatible with:**
 - M968FW Filter

2. Packaging:

- **Case:** One bag of 48 units

3. Product Specifications:

Compatible with the *Statspin Cytofuge-2 Cytocentrifuge*

- **Material:**
 - Polypropylene

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **LATEX:** Material is Latex Free
- **CE:** Product is CE marked

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination

8. Recommended Use:

- For use in **StatSpin Cytofuge 2 Cyto centrifuge** automated equipment, follow the equipment manufacturer's instructions for use

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*