

**Product Name:** Cryovial®

**Catalogue No.:** T301-5

**Material:**

Polypropylene Tube and Cap with Silicone seal

**Description:**

- 5 ml tube with cap
- Internal Threaded
- Round Bottom
- Printed with graduations and writing area
- Assembled with cap and red O-ring seal
- Compatible with color coding cap insert T312

**Packaging:**

**Bags:** 100 units

**Case:** 1000 units

**Standards and Conformity:**

**Gamma Sterilized:** SAL of  $10^{-3}$

**Radiation specified dose:** between 6.5 and 13.5 KGy

**Temperature range:** -196°C to +121°C

**Autoclavable:** +121°C for up to 30 minutes

**RNase, DNase, DNA and Pyrogen free:** certified

**Leak-proof:** tested in vacuum at 71.3cm Hg

**Liquid Nitrogen resistant:** Gas phase of LN<sup>2</sup>

**ISO 2859-1:** Sampling and inspection

**FDA:** Plastics conform to FDA 21 CFR 177.1520

**USP:** Plastic conforms to USP Class VI

**CONEG / RoHS:** Plastics comply

**REACH:** Plastic comply

**LATEX:** Plastic is Latex free

**BSE/TSE:** Plastic is BSE / TSE Free

**MEA:** Product has passed MEA testing

**Recommended Use:** Designed for storing biological material at temperatures as low as -196°C  
Not recommended in the liquid phase of LN<sup>2</sup>. Disposables for single use only.

**Storage Conditions:** Store at room temperature in normal warehouse conditions. Protect from any possible contamination and any damage to the packaging which could compromise the product sterility.

**Disclaimer:** *The information contained in this document shall not constitute a representation or warranty of the suitability for any specific purpose. To the best of our knowledge, the information contained is reliable and true based on our internal and external testing as well as statements supplied to us by third party suppliers.*