

**Product Name:** Paraffin Block Mailer

**Catalogue No.:** M477-6

This document replaces any previous version

**1. Product Description:**

- **Paraffin block mailer:** Non-sterile, 6 compartment container for transporting paraffin blocks.

**2. Packaging:**

- **Case:** 50 units per package

**3. Product Specifications:**

- **Material:**
  - PVC
- **Temperature range:** -30°C to +50°C
- **Autoclavable:** No
  
- **Compartment size:** 37 x 24 x 13.5 mm

**4. Standards and Conformity:**

- **ISO 2859-1:** Sampling and inspection procedures

**5. Quality Assurance:**

- No presence of contamination in plastic
- Visual attributes

**6. Traceability:**

- **Lot No. Composition:** 8 or 9 digits
  
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

**7. Storage Conditions:**

- Store at room temperature in normal warehouse conditions

**8. Recommended Use:**

- Accessory to be used in combination with cassettes for handling and shipping of paraffin blocks

## **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*