

Product Name: UniMailer™

Catalogue No.: M800-100 series

This document replaces any previous version

1. Product Description:

- **Slide mailer:** Exact fit to protect your slide during shipping. Fits one slide per tray. inter-locking base to stack, attach and ship multiple slides.
- **Available Colors:**
 - White, Pink, Green, Yellow and Blue
- **Compatible with:**
 - Fits standard 3 x 1 inch microscope slides

2. Packaging:

- **Case:** 4 packs of 50 units per case for a total of 200 units

3. Product Specifications:

- **Material:** High Impact Polystyrene
- **Temperature range:** -80°C to +80°C

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1640
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Assembly verification



6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging or product

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*