

Product Name: EasyDip™

Catalogue No.: M900-12 Series

This document replaces any previous version

1. Product Description:

- **Slide Staining Jar:** Non-sterile jar only
- **Available Colors:**
 - Blue, Green, Pink, White and Yellow. Also offered in assorted colors (1 of each color)
- **Compatible with:**
 - M905-12DGY Dark Gray Slide Staining Rack
 - M906 Aluminum Holder

2. Packaging:

- **Case:** 6 units per case (individually wrapped)

3. Product Specifications:

- **Material:** Acetal
- **Temperature range:** -80°C to +30°C
- **Can hold up to 12 slides with the staining rack**

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Avoid picric acid and phenol: prolonged use in these reagents will cause degradation
- Follow manufacturer's instructions when using colorant or other chemical substance

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*