

Product Name: **StainTray™**

Catalogue No.: **M918-1**

This document replaces any previous version

1. Product Description:

- **Slide staining system:** Slide staining and Humidity Chamber with raised rails to hold 10 slides
- **Compatible with:**
 - Fits standard 3 x 1 inch microscope slides

2. Packaging:

- **Case:** 1 individually wrapped unit per case

3. Product Specifications:

- **Material:**
 - Base: ABS plastic
 - Clear Cover: PETG
- **Temperature range:**
 - Base -80° C TO + 80° C
 - Cover -20° C TO + 60° C

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging or product

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*