

**Product Name:** **CytoSep™**

**Catalogue No.:** **M964-10FW1**

This document replaces any previous version

## 1. Product Description:

- **Single Cytology Funnel with Cap:** Disposable and pre-assembled with white paper filter

## 2. Packaging:

- **Case:** 5 packages of 100 (individually wrapped) units / 500 units per case

## 3. Product Specifications:

**Compatible with the Shandon Cytopspin® Cytocentrifuge:**

- **Material:**
  - Funnel: Polypropylene
  - Cap: Polypropylene
  - Filter: Cellulose Paper

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **LATEX:** Material is Latex Free
- **CE:** Product is CE marked

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag



## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination

## 8. Recommended Use:

- For use in **Shandon Cytospin® Cyto centrifuge** automated equipment, follow the equipment manufacturer's instructions for use

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*