

**Product Name: CytoSep™ Cytology Funnel****Catalogue No.: M965-20FW****This document replaces any previous version****1. Product Description:**

- **Cytology Funnel:** Disposable Double funnel pre-assembled with white paper filter and cap included

**2. Packaging:**

- **Case:** 12 inner boxes of 40 / case of 480

**3. Product Specifications:**

- Compatible with the Shandon Cytospin® Cytocentrifuge:
- Material:
  - Funnel: Polypropylene
  - Cap: Polyethylene
  - Filter: White Cellulose filter board

**4. Standards and Conformity:**

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Paper conforms to 21CFR 176.170 & 21CFR 176.180
- Plastic conforms to FDA 21 CFR 177.1520
- **CE:** Product is CE marked
- **LATEX:** Material is Latex Free
- **BSE / TSE:** Material is BSE / TSE Free

**5. Quality Assurance:**

- No presence of contamination in plastic or paper
- Assembly and fit
- Visual attributes
- Dimension measurements

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store in at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination

## 8. Recommended Use:

- For use in **Shandon Cytospin® Cytocentrifuge** automated equipment, follow the equipment manufacturer's instructions for use

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*