

Product Name: UniRack™ Jr.

Catalogue No.: S500-25 series

This document replaces any previous version

1. Product Description:

- **Rack:** Non-sterile, 25 position rack

- **Available colors:**
 - Blue, Red and Yellow

- **Compatible with:**
 - 1 to 5 ml Cryotubes and 0.5 to 2 ml microtubes
 - 0.2 to 0.5 ml PCR tubes

2. Packaging:

- **Case:** 10 units per case

3. Product Specifications:

- **Material:**
 - Polypropylene
- Temperature range: -196°C to +121°C
- Autoclavable at 121°C, for up to 30 minutes
- Total capacity: 25 cryovials and microtubes on one side or 16 PCR tubes on other side

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions;

8. Recommended Use:

Rack for holding tubes

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*