

**Product Name: Bioblock™
Deep Well Plate**

Catalogue No.: T110-200

This document replaces any previous version

1. Product Description:

- **200 ul capacity, deep well plate :** Non-sterile, 384 square wells with rounded bottom

Available colors:

Natural, *blue, *green, *pink and *yellow
*minimum quantity applies.

2. Packaging:

- **Case:** 4 packs of 6 units /24 units per case

3. Product Specifications:

- **Material:** Polypropylene
- **Certified RNase, DNase, Pyrogen and DNA Free**
- **Temperature range:** -196°C to +121°C
- **Autoclavable at +121°C for up to 30 minutes**
- **DMSO resistant**

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
- **CONEG / RoHS:** Plastic conforms to CONEG / RoHS standards for heavy metals
- **SBS standard:** Footprint conforms to SBS standard

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension and volume measurements.

6. Traceability:

- **Lot No. Composition :** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any damage to the packaging

8. Recommended Use:

- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*