

Product Name: Capinsert™

Catalogue No.: T345 Series

This document replaces any previous version

1. Product Description:

- **Color Coding Cap Insert:** Non-sterile, insert for color coding purposes
- **Available colors:** Blue, Gray, Green, Lilac, Orange, Pink, Red, Violet, White, Yellow
and Assorted (Blue, Lilac, Red, Yellow and White)
- **Compatible with:**
 - Microwtube® Cap Series
 - SecureTainer® and SpecTainer® Caps
 - T552-7WTP and T552-30WTP Caps

2. Packaging:

- **Case:** 500 units per bag

3. Product Specifications:

- **Material:**
 - Polypropylene
- Frosted Writing surface
- Temperature range: -196°C to +121°C
- Autoclavable at 121°C, for up to 30 minutes

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimensional measurements
- Insert verification

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions

8. Recommended Use:

- Should be used only in the gas phase of Liquid Nitrogen

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*