



## UNISLETTE™ 35° and 45° angle

# Tissue & Biopsy cassettes


Efficient flow-through slots maximize the fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. The snap latch and hinge-lock design of the UNISLETTE™ prevent early separation of base and lid and allows for one-hand operation. Resistant to the chemical action of most histological solvents.

### Specifications

- Anterior writing area at a 35° angle or **NEW** 45° angle
- Compatible with Leica printers
- Available in QuickLoad™ Stack (tape)
- Large front tab for easier lid opening
- Easy-to-remove disposable cover
- Biopsy cassettes have 1 mm square openings to maximize fluid exchange
- Made of acetal



# Ordering information



## For Leica printers

	35° degree angle			NEW 45° degree angle	
	Tissue	Biopsy	Color	Tissue	Biopsy
QuickLoad™ stacks 1000/Cs	M505-2T	M506-2T	○ White	M525-2T	M526-2T
	M505-3T	M506-3T	● Pink	M525-3T	M526-3T
	M505-4T	M506-4T	● Green	M525-4T	M526-4T
	M505-5T	M506-5T	● Yellow	M525-5T	M526-5T
	M505-6T	M506-6T	● Blue	M525-6T	M526-6T
	M505-7T	M506-7T	● Peach	M525-7T	M526-7T
	M505-8T	M506-8T	● Tan	M525-8T	M526-8T
	M505-9T	M506-9T	● Gray	M525-9T	M526-9T
	M505-10T	M506-10T	● Lilac	M525-10T	M526-10T
	M505-11T	M506-11T	● Orange	M525-11T	M526-11T
	M505-12T	M506-12T	● Aqua	M525-12T	M526-12T



Our original 35° angle cassette

Our NEW 45° angle cassette

Fluo Pink, Fluo Green and Fluo Yellow are also available. Minimum quantities apply. Contact us for more details.

## Accessories

### M474 Stainless Steel Base Molds

Manufactured from high-quality stainless steel for optimal thermal conductivity, the molds have a highly polished surface for easy paraffin block removal.

Cat. #	Size (mm)	Qty/Cs
M474-1	7 x 7 x 5	12
M474-2	15 x 15 x 5	12
M474-3	24 x 24 x 5	12
M474-4	30 x 24 x 5	12
M474-5	37 x 24 x 5	12



### M475 Disposable Base Molds

Inexpensive enough to be discarded after use, yet strong enough to be reused.

Cat. #	Size (mm)	Qty/Cs
M475-1	7 x 7 x 5	1000
M475-2	15 x 15 x 5	1000
M475-3	24 x 24 x 5	1000
M475-4	30 x 24 x 5	1000
M475-5	37 x 24 x 5	1000

