

| | |
|--|--|
| Product Name: Urine Collection Tube | Catalogue No.: T410-3 |
|--|--|

| |
|--|
| This document replaces any previous version |
|--|

1. Product Description:

- **12ml Urine Collection tube with Cap:** Disposable, non-sterile tube with cap.

2. Packaging:

- **Case:** 1 bag of 500 tubes and 5 bags of 100 caps per case

3. Product Specifications:

- **Material:**
 - Tube: Polystyrene
 - Lip Seal Cap: Low Density Polyethylene;
 - Cap Color: Blue
- Temperature range: -80°C to +80°C;
- Tubes have molded graduations from 1 ml to 12 ml;
- Centrifuge resistant up to 2 000 g.

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1640 and 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **BSE / TSE:** Material is BSE / TSE Free
- **LATEX:** Material is Latex Free

5. Quality Assurance:

- Clear, no presence of contamination in plastic;
- Visual attributes;
- Dimension and volume measurements;
- Closure verification.

6. Traceability:

- **Lot No. Composition:** 8 to 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label;
 2. On label inserted inside the master case;
 3. On the inner bag.

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions;
- Avoid temperature variations and humidity;
- Protect from any possible contamination;

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents;
- Follow chemical resistance chart recommendations;

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*