

Product Name: **Cultubes™**

Catalogue No.: **T425-6A**

This document replaces any previous version

1. Product Description:

- **5ml Culture Tube with Cap:** Disposable, sterile 12 x 75 mm culture tube, non-printed. Without Cap.

2. Packaging:

- **Case:** 20 bags of 125 units / 500 units per case

3. Product Specifications:

- **Material:** Polypropylene
- **Certified Pyrogen Free**
- **Gamma radiation sterilized** at a SAL of 10^{-3} ; specified dose between 6.5 kGy and 13.5 kGy
- **Temperature range tube:** -80°C to +121°C
- **Autoclavable**
- **Centrifuge resistant** at up to 3 000 g

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resins conform to FDA 21 CFR 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free
- **BSE / TSE:** Material is BSE / TSE Free

5. Quality Assurance:

- **Clear, no presence of contamination** in plastic
- **Visual attributes**
- **Dimension and volume measurements**
- **Closure verification**

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

➤ **The lot number can be found in one or all of these locations:**

1. On exterior case label
2. On label inserted inside the master case
3. On the inner bag

7. **Storage Conditions:**

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging which could compromise the product sterility

8. **Recommended Use:**

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*