

**Product Name: SlideFile™ Jr.**

**Catalogue No.: M700-50 series**

## 1. Product Description:

➤ **Slide storage system including removable rack:** Slides are stored in an upright position for easy insertion and removal. Rack has 50 numbered slots that fits 50 to 100 slides. For maximum storage space removal of the rack will allow to fit up to 200 slides.

➤ **Available Colors:**

- White, Pink, Green, Yellow and Blue

➤ **Compatible with:**

- Fits standard 3 x 1 inch microscope slides

## 2. Packaging:

➤ **Case:** 10 individually wrapped units per case

## 3. Product Specifications:

➤ **Material:**

- Base: High Impact Polystyrene
- Cover: Polystyrene

➤ **Temperature range:** -80°C to +80°C

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1640
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Closure verification

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
  
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging or product

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*