

**Product Name:** PCRAK™ with clear lid

**Catalogue No.:** T328-96 series

## 1. Product Description:

➤ **Rack:** Non-sterile, Stackable racks of 0.2 ml tubes assembled with transparent lid

➤ **Available colors:**

**Rack:** Blue, Green, Orange, Pink, Red Yellow and Assorted

**Lid:** Transparent

➤ **Compatible with:**

- 0.2 ml tubes, Strips of 8 or 12 tubes and 96-well PCR plates

## 2. Packaging:

➤ **Case:** 4 packs of 5 units per case

## 3. Product Specifications:

➤ **Material:**

- Polypropylene

➤ **Temperature range:** -80°C to +121°C

➤ **Autoclavable at 121°C,** for up to 30 minutes

➤ **Total capacity:** 96 wells

## 4. Standards and Conformity:

➤ **ISO 2859-1:** Sampling and inspection procedures

## 5. Quality Assurance:

➤ Clear, no presence of contamination in plastic

➤ Visual attributes

➤ Dimension measurements

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
  
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions;

## 8. Recommended Use:

Rack for holding 2.0 ml tubes, Strips of 8 or 12 tubes and 96-wells PCR plates

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*