

TECHNICAL DATA SHEET

EDL41218

PRODUCT CODE

EDL41218

PACKAGING

Cases of 250

SHELF LIFE

Five years sterility guaranteed

PRODUCT DESCRIPTION

Sterile sampling bag made of FDA approved virgin polyethylene tubing with tear off and easy pull tabs to open the bag. Bags are easily closed thanks to two flat wires covered with white strips. Product is folded for convenient storage.

PRODUCT SPECIFICATIONS

Closure type	2 FLAT WIRES with Safety Tabs
Closure material	Steel wires covered with FDA approved pressure sensitive tape with a backing of polypropylene or PVC.
Printed	Clear – non printed material
MSDS	Not applicable for this product
Sterilisation	Sterility is assured due to extrusion temperatures exceeding 220°C which eliminates any present bacteria. During manufacturing all bags are sealed at both top and bottom in a continuous format before being cut into individual bags. This method assures that the bag's interior will never come in contact with its surroundings. Periodic validation tests, based on USP 40 – Sterility testing, are performed in our facility to guaranty the sterility of our processes.
RNase/DNase/Pyrogens	Products are free of RNase, DNase and pyrogens
Volume	228 oz (6743 ml)
Dimensions	12 X 18 in (305 mm x 457 mm)
Thickness of material	4 mil (102 micron)
FDA regulation	All TWIRL'EM sterile sampling bags comply to the applicable regulations of the Food and Drug Administration (21 CFR 177.1520 (c) 3.1a, (c)3.2a and 178.2010.
EU regulation	All TWIRL'EM sterile sampling bags comply where applicable to EU No 10/2011 regulation for use in contact with food, subject to certain types of migration limitations and labeling requirements.
CFIA	All TWIRL'EM sterile sampling bags are also a solution which could be used in the CFIA Prevention Control Plan inspired by HACCP principles, a Canadian federal initiative from the Safe Food for Canadians Regulations (SFCR).