

Technical Data Sheet

Product Name: Cryosette™

Catalogue No.: M9560

This document replaces any previous version

1. Product Description:

- **Tissue Storage Container with Lid:** Disposable and non-sterile
- **Available Colors:**
 - Orange container
 - White lid

2. Packaging:

- **Case:** 250 containers with 5 bags of 50 lids per case / 250 units per case

3. Product Specifications:

- **Material:**
 - Acetal container
 - High Density Polyethylene lid
- Slanted Writing surface
- Temperature range: -196°C to +30°C

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.2470 and CFR 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- closure verification

Technical Data Sheet

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging

8. Recommended Use:

- Verify proper lid closure when using biohazard material
- Avoid picric acid, prolonged use in this acid will cause degradation
- Follow chemical resistance chart recommendations

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*