


**M970** *Series*

*Simport*

*Eliminate transcription error.  
Complete only 1 identification label instead of 12!*

# The **CoreDish**<sup>®</sup> Multiple Biopsy Containers

Requires only 30ml of Formaldehyde instead of 120ml.  
12 specimens in **1 leakproof container.**



**Only one leakproof closure to open and close, not 12.**

A fail proof concept, placing the biopsy in the proper compartment to guarantee no mix up or error.  
This design prevents any tissue from migrating from one compartment to another.  
Almost 50% less plastic than 12 containers.

If you **REALLY** care about your biopsies, let us help you **FACILITATE** their handling!


Prefiled with 10% formalin or without




[simport.com](http://simport.com)

## Protect Tissue Integrity: All-in-one system eliminating excessive handling of delicate biopsy specimens.

The Simport® CoreDish® is a multi-compartment leakproof container in the shape of a dish for holding and transporting biopsies. Its closure is equipped with a silicone o-ring ensuring total protection of contents. The Simport® CoreDish® measures only 15 x 95 mm in diameter. Each compartment is clearly identified to allow proper placement and visualization of the biopsy being inserted. Thanks to the CoreDish® it is no longer necessary to identify and use a multitude of individual containers, thereby reducing the risks of confusion during their manipulation or identification of patient information. The Simport® CoreDish® offers many configurations in order to hold different biopsies of the breast, prostate, upper GI tract and lower GI tract. A unique label allows for all essential information to be written such as patient I.D., doctor, date and time. To ensure complete integrity of contents from the collection stage to the reopening of the CoreDish® a security label that seals the unit is supplied.



Cat. #	For	Compartments	Prefilled	Qty/Pk	Qty/Cs
M970-D5B-1	Breast	5	NO	1	10
M970-D5B-2	Breast	5	YES	1	10
M971-D5B-2	Breast	5	NO	1	10
M970-D8P	Prostate	8	YES	1	10
M971-D8P	Prostate	8	NO	1	10
M970-D8UGI	Upper Gi Tract	8	YES	1	10
M971-D8UGI	Upper Gi Tract	8	NO	1	10
M970-D12LGI	Lower Gi Tract	12	YES	1	10
M971-D12LGI	Lower Gi Tract	12	NO	1	10
M970-D12P	Prostate	12	YES	1	10
M971-D12P	Prostate	12	NO	1	10



A security label is supplied to ensure integrity of contents from collection stage to reopening of CoreDish®.

To close and seal, turn clockwise until you feel a firm stop.

### M975 The CorePicker™ Made of polystyrene

A practical tool and a great help to pick up and handle biopsies out of the CoreDish®. Packed in tamperproof resealable bags.

Cat. #	Length	Qty/Bag	Qty/Pk
M975	53 mm (2 1/8 in.)	25	125

### M976 Shipping Box for CoreDish® Made of cardboard

This sturdy and easy-to-assemble shipping box is most convenient for transporting or mailing the Simport® CoreDish®.

Cat. #	Description	Qty/Pk
M976	Cardboard Box	10