

# Technical Data Sheet

**Product Name:** SuperFLO

**Catalogue No.:** M535

## Micro Biopsy Cassette with assembled hinged Lid

**Material:** Acetal

**Case:** 500 units per pack, 1500 per case

**Description:** Disposable, Non-sterile  
With assembled lid, a 45° angle and a surface ideal for handwriting

**Temperature range:** -80°C to +30°C

**ISO 2859-1:** Sampling and inspection

**FDA:** Plastics conform to FDA 21 CFR 177.2470

**CONEG / RoHS:** Plastics and colorants comply

**LATEX:** Material is Latex free

**CE:** Product is CE marked

**Recommended Use:** For histological sample containment and processing

**Storage Conditions:** Store at room temperature in normal warehouse conditions. Protect from any possible contamination and any damage to the packaging.

**OWNERSHIP OF MATERIALS:** *Materials and information contained in this document are Simport Scientific's copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*