

# Technical Data Sheet

**Product Name:** SuperFLO

**Catalogue No.:** M535-SL

## Micro Biopsy Cassette with separate hinged Lid

**Material:** Acetal

**Case:** 10 sleeves of 75 cassettes with 10 bags of 75 lids / 750 Total units

**Description:** Disposable, Non-sterile  
With separate lid, a 45° angle and compatible with Thermo Fisher cassette printers

**Temperature range:** -80°C to +30°C

**ISO 2859-1:** Sampling and inspection

**FDA:** Plastics conform to FDA 21 CFR 177.2470

**CONEG / RoHS:** Plastics and colorants comply

**LATEX:** Material is Latex free

**CE:** Product is CE marked

**Recommended Use:** For histological sample containment and processing

**Storage Conditions:** Store at room temperature in normal warehouse conditions. Protect from any possible contamination and any damage to the packaging.

**OWNERSHIP OF MATERIALS:** *Materials and information contained in this document are Simport Scientific's copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*