

<b>Product Name: CRYO CUBE STORAGE BOX</b>	<b>Catalogue No.: T315-81Y</b>
--	--------------------------------

## CRYO CUBE STORAGE BOX YELLOW

**Material:** Box: Polypropylene

**Case:** 8 units per case

**Description:** Reusable Polypropylene (PP) Cryo Cube Box/ Storage Cube Box  
Square Box with Transparent Lid.

**Temperature range:** -80°C to +121°C

**Pass ISO 10993:** External communication devices for indirect blood contact for prolonged period.

**Pass USP Class IV**

**FDA:** Plastics conform to US FDA 21 CFR

**CONEG / RoHS:** Plastics comply

**BSE/TSE free**

**LATEX:** Material is Latex free

**Recommended Use:** Box for storage holding up to 81 tubes of up to 2mL capacity.

**Storage Conditions:** Store at room temperature in normal warehouse conditions. Protect from any possible contamination and any damage to the packaging.

**OWNERSHIP OF MATERIALS:** *Materials and information contained in this document are Simport Scientific's copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*